

**Adjustable  
Modular  
Spinal System**



# **SpinaLoc** Spinal Orthosis

The Spinal Technology SpinaLoc offers unsurpassed fit, functionality, and flexibility in a prefabricated spinal brace. It provides firm control of the torso, eliminating unwanted movement to stabilize the spine, as well as a greater degree of patient comfort and pain abatement.

The SpinaLoc achieves both patient comfort and firm control by integrating rigid elements and a soft, flexible patient interface. The modular system makes fitting procedures fast and effective. The SpinaLoc provides anterior, posterior and lateral control. The anterior and posterior components can be further adjusted by simply heating and contouring the inserted panels.

For more information and pricing, contact your regional Sales Representative or call Customer Service at:  
**1 800 253 7868.**

## **Features**

**Non-allergenic foam allows for pressure relief and helps reduce skin breakdown**

**Rigid plastic panel inserts provide firm control**

**A patient friendly dynamic closure system with extensive adjustability**

## **SpinaLoc** Standard Orthosis **Anterior View**



## **Posterior View** Standard Orthosis

## **Features**

**Removable Posterior and Anterior panels that can be modified to patients' diagnoses**

**Removable Lateral support panel for stability**



## **SpinaLoc II** Low Profile Orthosis

### **Anterior View**



### **Posterior View**



**Spinal**  
Technology, Inc.





# SpinaLoc

## Spinal Orthosis

### Adjustable Modular Spinal System

- A patient-friendly dynamic closure system with extensive adjustability
- Removable lateral support panel for stability
- Non-allergenic foam allows for pressure relief and helps reduce skin breakdown
- Removable anterior and posterior panels that can be modified to patients diagnoses
- Sternal shield to control thoracic area
- Thoracic anti-kyphotic support
- Full abdominal panels

Size Chart

### SpinaLoc

| Men's<br>Size | Hip, Waist &<br>Xyphoid<br>Circumference | Anterior<br>Height | Posterior<br>Height |
|---------------|--|--------------------|---------------------|
| M1 (S)        | 26" - 40"                                | 11"                | 16"                 |
| M2 (M)        | 32" - 46"                                | 12"                | 16"                 |
| M3 (L)        | 40" - 54"                                | 13"                | 16"                 |
| M4 (XL)       | 44" - 58"                                | 13"                | 18"                 |

### Women's

| Size      |           |       |     |
|-----------|-----------|-------|-----|
| F1 (S/M)  | 25" - 37" | 9.5"  | 15" |
| F2 (M/L)  | 33" - 45" | 10.5" | 15" |
| F3 (L/XL) | 41" - 53" | 11.5" | 16" |

Low Profile

|     |           |      |     |
|-----|-----------|------|-----|
| F1L | 25" - 37" | 6.5" | 15" |
| F2L | 33" - 45" | 7"   | 15" |
| F3L | 41" - 53" | 7.5" | 16" |

### SpinaLoc II

| Size     |           |      |       |
|----------|-----------|------|-------|
| Small    | 26" - 34" | 6.5" | 10.5" |
| Medium   | 31" - 39" | 6.5" | 11"   |
| Large    | 36" - 44" | 7"   | 11"   |
| Ex-Large | 41" - 49" | 7"   | 12"   |
| XX-Large | 46" - 55" | 7.5" | 12"   |
| 3X-Large | 52" - 60" | 7.5" | 12"   |

### Additional SpinaLoc Components Available

#### Sternal Shield

##### Anterior View

with Anterior Thoracic Extension (ATE)

#### Features

Optional  
Sternal Shield  
available  
with easy  
attachment and  
adjustable  
straps



with Shoulder Straps

with Axilla Straps

#### Thoracic Extension

##### Posterior & Anterior View

with Posterior Thoracic Extension (PTE)

Thoracic  
Extension &  
Anti-Kyphotic  
component  
to control  
Thoracic region



#### Inserts

Additional  
Inserts  
for lateral  
stability are  
also  
available



Spinal  
Technology, Inc.

