

LOOK WHAT'S NEW

Spinal
Technology, Inc.



**Prefabricated
Spinal System**

S.T.O.P. Orthosis
Spinal Trauma Orthotic Positioning

SADMERC Classified

The S.T.O.P. orthosis is an extremely versatile spinal immobilization device fabricated in multiple configurations.

- Bi-valve with molded tongues for maximum adjustability
- Anatomically contoured to facilitate superior fit on a wide range of patients
- Removable washable lining for increased comfort and compliance
- Available rigid molded side panels to increase lateral stability

S.T.O.P. I & II (2 piece design)

- S.T.O.P II is MRI compatible
- Female design with breast cut out



L0460

L0462



S.T.O.P. III (3 piece design)

- Adjustable Anterior Thoracic Extension (ATE)

S.T.O.P. IV (4 piece design)

- Adjustable Anterior Thoracic Extension (ATE)
- Adjustable Posterior Thoracic Extension (PTE)



L0464

Spinal
Technology, Inc.



800 253 7868
888 775 0588 / Fax

S.T.O.P. Orthosis

Spinal Trauma Orthotic Positioning

Spinal
Technology, Inc.



Prefabricated
Spinal System

Size Chart

Male Sizes	Hip	Waist	Chest
Neutral Abdominal			
Small	35" - 38"	32" - 35"	35" - 38"
Medium	39" - 42"	36" - 39"	39" - 42"
Large	43" - 46"	40" - 43"	43" - 46"
X-Large	47" - 50"	44" - 47"	47" - 50"
Abdominal Relief			
Small	35" - 38"	34" - 37"	35" - 38"
Medium	39" - 42"	38" - 41"	39" - 42"
Large	43" - 46"	42" - 45"	43" - 46"
X-Large	47" - 50"	46" - 49"	47" - 50"
Female Sizes			
Neutral Abdominal			
Small	33" - 36"	29" - 32"	31" - 34"
Medium	37" - 40"	33" - 36"	35" - 38"
Large	41" - 44"	37" - 40"	39" - 42"
X-Large	45" - 48"	41" - 44"	43" - 46"
Abdominal Relief			
Small	33" - 36"	31" - 34"	31" - 34"
Medium	37" - 40"	35" - 38"	35" - 38"
Large	41" - 44"	39" - 42"	39" - 42"
X-Large	45" - 48"	43" - 46"	43" - 46"

Spinal
Technology, Inc.



800 253 7868
888 775 0588 / Fax